Tri Auto ZX2+

Thinking ahead. Focused on life.

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Motore Endodontico Intelligente

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* Decide on the final file size according to the shape of the root canal.

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

Casi

Case 1

Curved root canal (left maxillary first molar)


Preoperative

Postoperative

Mesiobuccal (MB) and distobuccal (DB) roots of #26 are constricted, and only a few millimeters of patency was achieved with hand files. Then a mechanical glide path using OGP2 was established. OGP2 could obtain glide path preparation without any difficulties such as broken files or ledge formation to constricted canals even to DB, which has severe constriction and S-shape curvature. Total time to complete mechanical glide path using OGP2 was 5 minutes.

case 2

Sclerosed root canal (right maxillary first molar)

Preoperative

Postoperative

The radiographic images highlight that a root canal system was still present. The detection of the root canal orifice and root canal preparation were difficult due to the presence of mineralized tissue. The OGP2 mode contributed to attaining patency and glide path on 4 roots, including the second mesiobuccal (MB2) root canal. Subsequently, root canal preparation was completed using OTR mode on the device. All procedures ranging from the patency of the 4 root canals to the final root preparation were carried out in approximately 7 minutes.

case 3

Patency abandoned (left maxillary second molar)

Preoperative

Postoperative

Retreatment case of #37, root canal patency was not obtained during the past treatment. After root canal orifices are specified, patency was obtained using OGP2. Both mesiobuccal (MB) and distobuccal (DB) canals have S-shape curvature, and especially MB was challenging. If only with hand instruments, treatment time for MB would be extremely long. By using OGP2, total time needed for glide path preparation was 14 minutes in total; 6 minutes for MB, 5 minutes for DB and 3 minutes for palatal canal.

Clinical images and comments are provided by: Dr. Toshihiro Ushikubo
Dr. Motoki Okamoto
Dr. Tai Gega

*The images shown here represent examples of post surgery progress at the clinic that provided the photos.

*A pre-sales clinical survey conducted by J.MORITA MFG.CORP. showed that of the dentists that used OGP2, 77% used it for patency and 86% used it for glide path.

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Accessories: Tester, LS OIL
 Power Plugs (one each of 4 types)
 Probe Cord (0.75 m)
 File Holder, Contrary Electrodes
 HP Protective Sleeve Type A

Options: Handpiece Holder, Long File Holder
 External File Electrode (with cap)
 Probe Cord (1.8 m)

Handpiece


Free Running Operation
 Speed: 100 – 1000 rpm
 Torque: 4 Ncm or more
 Battery: Lithium-ion battery (DC 3.7 V)
 Weight: Approx. 140 g (including battery and contra angle)

Battery Charger

Rated Input Voltage: DC 5 V
 Rated Input Current: 2.4 A
 Weight: Approx. 280 g
 (AC adapter not included)

AC Adapter

Rated Input Voltage: AC100 - 240 V
 Rated Input Frequency: 47 - 63 Hz
 Rated Input Current: 0.4 A

Development and Manufacturing

J. MORITA MFG. CORP.

680 Higashihama Minami-cho, Fushimi-ku,
Kyoto 612-8533, Japan
T +81. (0)75. 611 2141, F +81. (0)75. 622 4595

Morita Global Website

www.morita.com

Distribution

J. MORITA CORP.

3-33-18 Tarumi-cho, Suita-shi, Osaka 564-8650, Japan
T +81. (0)6. 6380 1521, F +81. (0)6. 6380 0585

J. MORITA USA, INC.

9 Mason, Irvine CA 92618, USA
T +1. 949. 581 9600, F +1. 949. 581 8811

J. MORITA EUROPE GMBH

Justus-von-Liebig-Strasse 27b, 63128 Dietzenbach, Germany
T +49. (0)6074. 836 0, F +49. (0)6074. 836 299

MORITA DENTAL ASIA PTE. LTD.

150 Kampong Ampat, #06-01A KA Centre, Singapore 368324
T +65. 6779. 4795, F +65. 6777. 2279

J. MORITA CORP. AUSTRALIA & NEW ZEALAND

Suite 2.05, 247 Coward Street, Mascot NSW 2020, Australia
T +61. (0)2. 9667 3555, F +61. (0)2. 9667 3577

J. MORITA CORP. MIDDLE EAST

4 Tag Al Roasaa, Apartment 902, Saba Pacha 21311 Alexandria, Egypt
T +20. (0)3. 58 222 94, F +20. (0)3. 58 222 96

J. MORITA CORP. INDIA

Filix Office No.908, L.B.S. Marg, Opp. Asian Paints, Bhandup (West), Mumbai 400078, India
T +91-22-2595-3482

J. MORITA MFG. CORP. INDONESIA

28F, DBS Bank Tower, Jl. Prof. Dr. Satrio Kav. 3-5, Jakarta 12940, Indonesia
T +62-21-2988-8332, F + 62-21-2988-8201

SIAMDENT CO., LTD.

71/10 Mu 5, Thakham, Bangpakong, Chachuengsao 24130, Thailand
T +66. 38. 573042, F +66. 38. 573043
www.siamdent.com

Diagnostic and Imaging Equipment

Treatment Units

Handpieces and Instruments

Endodontic System

Laser Equipment

Laboratory Devices

Educational and Training Systems

Auxiliaries

Morita Italia

via Villa 17\b Bresso Mi
T: 02.87394532 www.jmoritaitalia.com